# Automated Functional Imaging (AFI)

Leading the way in strain imaging  
from research to clinical routine

White Paper

**Uwe Lempetz**

Global Clinical Development Director

**Gunnar Hansen**

Global Clinical Research Manager

**Sten Roar Snare**

Engineering Leader Display & Image Analysis

The word 'vivid' is written in a large, white, sans-serif font. Behind the letters, there is a colorful, wavy graphic element that transitions through a spectrum of colors including green, yellow, orange, red, and purple, creating a sense of motion and energy.

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# References:

- 1 Zito et.al. Ten Years of 2D Longitudinal Strain for Early Myocardial Dysfunction Detection: A Clinical Overview. BioMed Research International, Volume 2018, Article ID 8979407
- 2 Plana et.al. Expert Consensus for Multimodality Imaging Evaluation of Adult Patients during and after Cancer Therapy: A Report from the American Society of Echocardiography and the European Association of Cardiovascular Imaging. J Am Soc Echocardiogr. 2014 Sep;27(9):911-39
- 3 Oxford Analytica; Automated Function Imaging Healthymagination Validation Summary. Oxford Analytica 2012
- 4 Reant et.al, Evaluation of Global Left Ventricular Systolic Function Using Three-Dimensional Echocardiography Speckle-Tracking Strain Parameters. J Am Soc Echocardiogr. 2012 Jan;25(1):68-79.
- 5 Li et.al. Value of Three-Dimensional Speckle-Tracking in Detecting Left Ventricular Dysfunction in Patients with Aortic Valvular Diseases. J Am Soc Echocardiogr. 2013 Nov;26(11):1245-52.
- 6 Uusitalo et al, Two-Dimensional Speckle-Tracking during Dobutamine Stress Echocardiography in the Detection of Myocardial Ischemia in Patients with Suspected Coronary Artery Disease. J Am Soc Echocardiogr 2016;29:470-9.
- 7 Brunet-Bernard et.al. Defining patients at-risk of non-response to cardiac resynchronization therapy. Value of rest and exercise echocardiography. International Journal of Cardiology, doi: 10.1016/j.ijcard.2013.11.085
- 8 Wang et al, Left Ventricular Long-Axis Performance during Exercise is An Important Prognosticator in Patients with Heart Failure and Preserved Ejection Fraction. Int J Cardiol 2015 Jan 15;178:131-5.
- 9 Galli et al, Value of Myocardial Work Estimation in the Prediction of Response to Cardiac Resynchronization Therapy. J Am Soc Echocardiogr. 2018 Feb;31(2):220-230.
- 10 Edwards et al, Global Myocardial Work Is Superior to Global Longitudinal Strain to Predict Significant Coronary Artery Disease in Patients With Normal Left Ventricular Function and Wall Motion. J Am Soc Echocardiogr. 2019 Aug;32(8):947-957.
- 11 El Mahdiui et al, Global Left Ventricular Myocardial Work Efficiency in Healthy Individuals and Patients with Cardiovascular Disease. J Am Soc Echocardiogr. 2019 Sep;32(9):1120-1127.
- 12 Pihadi et al, Determinants and prognostic implications of left ventricular mechanical dispersion in aortic stenosis. European Heart Journal - Cardiovascular Imaging 2019 Jul 1;20(7):740-748.
- 13 Badano et al, Standardization of left atrial, right ventricular, and right atrial deformation imaging using two-dimensional speckle tracking echocardiography: a consensus document of the EACVI/ASE/ Industry Task Force to standardize deformation imaging Eur Heart J -Cardiovasc Imaging 2018;19, 591-600
- 14 Thomas et al, Left Atrial Structure and Function, and Left Ventricular Diastolic Dysfunction. J Am Coll Cardiol 2019;73:1961-77;
- 15 Florescu et al, Automated left atrial volume measurement by two-dimensional speckle-tracking echocardiography: feasibility, accuracy, and reproducibility. European Heart Journal - Cardiovascular Imaging 2021; 00, 1-10
- 16 Sanz et al, Anatomy, Function, and Dysfunction of the Right Ventricle. Journal of the American College of Cardiology, Volume 73, Issue 12, 2 April 2019, Pages 1463-1482
- 17 Jones et al., Echocardiographic Assessment of the Right Ventricle—State of the Art. Heart,Lung and Circulation, Volume 28, Issue 9, 1339 – 1350
- 18 Voigt et al, Definitions for a Common Standard for 2D Speckle Tracking Echocardiography: Consensus Document of the EACVI/ASE/Industry Task Force to Standardize Deformation Imaging. J Am Soc Echocardiogr 2015;28:183-93
- 19 D’Hooge et al, Two-dimensional speckle tracking echocardiography: standardization efforts based on synthetic ultrasound data. Eur Heart J Cardiovasc Imaging 2016 Jun;17(6):693-701
- 20 Farsalinos et al., Head-to-Head Comparison of Global Longitudinal Strain Measurements among Nine Different Vendors The EACVI/ASE Inter-Vendor Comparison Study J Am Soc Echocardiogr 2015;28:1171-81
- 21 Mirea et al, Intervendor Differences in the Accuracy of Detecting Regional Functional Abnormalities. Eur Heart J Cardiovasc Imaging 2016 Jun;17(6):693-701.
- 22 Ünlü et al, Comparison of Feasibility, Accuracy, and Reproducibility of Layer-Specific Global Longitudinal Strain Measurements Among Five Different Vendors: A Report from the EACVI-ASE Strain Standardization Task Force. J Am Soc Echocardiogr 2019;32:624-32
- 23 Ünlü et al, Impact of apical foreshortening on deformation measurements: a report from the EACVI-ASE Strain Standardization Task Force European Heart Journal -Cardiovascular Imaging 2019; 0, 1-7

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